



Circumferential resection accessories



Selectable resection accessories



C-arm probe allows hemorrhoids treatment for bigger piles while prevents severe stricture issues from PPH.

Application:

1. Titanium alloy 3D staples provide superior strength to prevent anastomotic fistula.
2. A slim designed anvil puts less strain on the anastomotic site during removal.
3. Two staggered row of staples allows a good anastomosis and haemostatic seal.
4. Firing in the green range delivers a closed staple height appropriate for compressed tissue thickness.
5. Extra-large cartridge volume for more tissue incorporation.
6. Variety of sizes provide more options for the individual patient anatomy requirements.
7. Non-Slip Grip provides improved comfort and control.
8. Audible and visible feedback for better reminding surgeons about fire completion.

Model	Cartridge Color	Staple height (mm)	Closed staple height (mm)	Staple quantity (pc)	Anvil diameter (mm)
PPH33	(Green)	4.0	0.5-1.0	32	33.4
PPH34	(Blue)	4.0	0.5-1.0	32	33.6
PPH36	(White)	4.0	0.5-1.0	34	35.7

We offer more options for disposable anorectal staplers. Please consult to sales department for further information.